

## Retroviral gene transfer protocol

Barton lab

7/24/07

### Day 0:

Split GP2-293T/ $\Phi$ x packaging cell line for transfection the following day:

6cm plate:  
1.75-2.0x10<sup>6</sup> cells

10cm plate:  
7.5x10<sup>6</sup> cells

### Day 1:

Transfect packaging cell line by CaPO<sub>4</sub>

6cm plate:  
5 $\mu$ g DNA  
3 $\mu$ g vsv-g

10cm plate:  
15-20 $\mu$ g DNA  
10 $\mu$ g vsv-g

For CaPO<sub>4</sub> transfection:

**Pre-incubate** packaging cells 30' with **25 $\mu$ M chloroquine**

6cm plate:  
1.5 $\mu$ L 50mM stock  $\rightarrow$  3mL

10cm plate:  
5 $\mu$ L 50mM stock  $\rightarrow$  10mL

for ea plate combine *in numbered order* in a 15mL conical:

**1)** 8-16 $\mu$ g DNA

25-50 $\mu$ g DNA

**2)** H<sub>2</sub>O to 0.438 mL

H<sub>2</sub>O to 1.314 mL

**3)** 62 $\mu$ L CaCl<sub>2</sub>

186 $\mu$ L CaCl<sub>2</sub>

using a 10mL pipette bubble DNA and add drop-wise (as slow as you can stand):

0.5mL 2x HBS

1.5mL 2x HBS

**replace media** after ~8hrs:

remove 2mL, add 3mL

remove 8mL, add 10mL

incubate o/n @ 37°

### Day 2:

Replace media on transfection

6cm plate:  
remove 3mL, add 2mL

10cm plate:  
remove 10mL, add 6mL

**transfer** to 32° incubator o/n

split target cell line for transduction the next day:

6 well plate:  
**293T** 3x10<sup>5</sup> cells  
**3T3** 1.0x10<sup>5</sup> cells  
**RAW** 3x10<sup>5</sup> cells

10cm plate:  
18.0x10<sup>5</sup> cells  
6.0x10<sup>5</sup> cells  
18.0x10<sup>5</sup> cells

**Day 3:**

Harvest virus supernatant:

Filter sups through a 10cc syringe fitted w/ 22-45 $\mu$ m filter into a 15 mL conical

*Optional virus concentration step from 10cm plate:*

*Transfer sup to tubes for Sw40 Ti swinging bucket rotor*

*Spin at 50,000g (16.8rpm) for 90' in ultracentrifuge*

*Suck sup off of viral pellet*

*Re-suspend in 1mL media (1mL viral sup  $\rightarrow$  1 well of target cells in a 6 well plate in 2mL media)*

**Add polybrene** (hexadimethrine bromide) @ 5 $\mu$ g/mL to viral sup  
(1 $\mu$ L/mL of 5mg/mL stock)

Suck supernatant off of target line, apply retroviral sup

Incubate o/n at 32°C

**Day 4:**

Suck off viral sup and replace media on the infected cells and return to 37°C

**Day 5:**

Analyze cells for expression of marker